

Next generation plasma collection technology for the clinical analysis of temozolomide by HILIC/MS/MS

ASMS 2014 WP641

Alan J. Barnes¹, Carrie-Anne Mellor²,
Adam McMahon², Neil Loftus¹

¹Shimadzu, Manchester, UK

²WMIC, University of Manchester, UK

Next generation plasma collection technology for the clinical analysis of temozolomide by HILIC/MS/MS

Introduction

Plasma extraction technology is a novel technique achieved by applying a blood sample to a laminated membrane stack which allows plasma to flow through the asymmetric filter whilst retaining the cellular components of the blood sample.


Plasma separation card technology was applied to the quantitative analysis of temozolomide (TMZ); an oral imidazotetrazine alkylating agent used for the treatment of Grade IV astrocytoma, an aggressive form of brain tumour.

Under physiological conditions TMZ is rapidly converted to 5-(3-methyltriazene-1-yl)imidazole-4-carboxamide (MTIC) which in-turn degrades by hydrolysis to 5-aminoimidazole-4-carboxamide (AIC). Storage of plasma has previously shown that both at -70C and 4C degradation still occurs. In these experiments, whole blood containing TMZ standard was applied to NoviPlex plasma separation cards (PSC). The aim was to develop a robust LC/MS/MS quantitative method for TMZ.

Materials and Methods

Plasma separation

TMZ spiked human blood calibration standards (50uL) were applied to the PSC as described below in figure 1.



A NoviPlex card is removed from foil packaging.

Approximately 50uL of whole blood is added to the test area.

After 3 minutes, the top layer is completely removed (peeled back).

The collection disc contains 2.5uL of plasma. Card is air dried for 15 minutes.

Figure 1. NoviPlex plasma separation card workflow

Next generation plasma collection technology for the clinical analysis of temozolomide by HILIC/MS/MS




Figure 1. Noviplex plasma separation card workflow (Cont'd)




Figure 2. Applying a blood sample, either as a finger prick or by accurately measuring the blood volume, to the laminated membrane stack retains red cells and allows a plasma sample to be collected. The red cells are retained by a combination of adsorption and filtration whilst plasma advances through the membrane stack by capillary action. After approximately three minutes the plasma Collection Disc was saturated with an aliquot of plasma and was ready for LC/MS/MS analysis.

Sample preparation

TMZ was extracted from the dried plasma collection discs by adding 40uL acetonitrile + 0.1% formic acid, followed by centrifugation 16,000g for 5 min. 30uL supernatant was added directly to the LC/MS/MS sample vial for analysis.

As a control, conventional plasma samples were prepared by centrifuging the human blood calibration standards at 1000g for 10min. TMZ was extracted from 2.5uL of plasma using the same extraction protocol as applied for PSC.

Next generation plasma collection technology for the clinical analysis of temozolomide by HILIC/MS/MS

LC/MS/MS analysis

Ionisation	: Electrospray, positive mode MRM 195.05 >138.05 CE -10
Desolvation line	: 300°C
Drying/Nebulising gas	: 10L/min, 2L/min
Heating block	: 400°C
HPLC	: HILIC Nexera UHPLC system
Flow rate	: 0.5mL/min (0-7min), 1.8mL/min (7.5min-17.5min)
Mobile phase	: A= Water + 0.1% formic acid B= Acetonitrile + 0.1% formic acid
Gradient	: 95% B – 30%% B (6.5 min), 30% B (7.5 min), 95% B (18 min)
Analytical column	: ZIC HILIC 150 x 4.6mm 5um 200 ^a
Column temperature	: 40°C
Injection volume	: 10uL
	Reverse Phase Nexera UHPLC system
	0.4mL/min
	A= Water + 0.1% formic acid B= methanol + 0.1% formic acid
	5% B – 100%% B (3 min), 100% B (7 min), 5% B (10 min)
	Phenomenex Kinetex XB C18 100 x 2.1mm 1.7um 100A
	50°C
	2µL

Results

HILIC LC/MS/MS

Temozolomide is known to be unstable under physiological conditions and is converted to 5-(3-methyltriazen-1-yl)imidazole-4-carboxamide (MTIC) by

a nonenzymatic, chemical degradation process. Previous studies have used HILIC to analyze the polar compound and to avoid degradation in aqueous solutions.




Figure 3. HILIC LC/MS/MS chromatograms for PSC TMZ analysis at 0.5 and 8ug/mL. The PSC calibration curve was linear between 0.2-10ug/mL ($r^2 > 0.99$). HILIC was considered in response to previous published data and to minimize potential stability issues. However, to reduce sample cycle times a reverse phase method was also developed.

Next generation plasma collection technology for the clinical analysis of temozolomide by HILIC/MS/MS

Reversed Phase LC/MS/MS




Figure 4. Reverse phase LC/MS/MS chromatograms for PSC TMZ analysis at 0.5 and 8ug/mL. The PSC calibration curve was linear between 0.2-10ug/mL ($r^2 > 0.99$; replicate samples were included in the liner regression analysis at 0.5 and 8ug/mL; n=3).

Due to the relatively long cycle time (18 min), a faster reversed phase method was developed (10 min). Sample preparation was modified with PSC sample disk placed in 40uL methanol + 0.1% formic acid, followed by centrifugation 16,000g 5 min. 20uL supernatant was

added directly to the LC/MS sample vial plus 80uL water + 0.1% formic acid. In addition to reversed phase being faster, the sample injection volume was reduced to just 2uL as a result of higher sensitivity from narrower peak width (reversed phase, 13 sec; HILIC, 42 sec).

Comparison between PSC and plasma




Figure 5. Human blood TMZ calibration standards were prepared using PSC and conventional plasma. Using the confirmatory ion transition 195.05>67.05 both the PSC and plasma sample are in broad agreement with regard to matrix ion signal distribution.

Next generation plasma collection technology for the clinical analysis of temozolomide by HILIC/MS/MS




Figure 6. Human blood TMZ calibration standards were prepared using PSC and conventional plasma. Using the quantitation ion transition 195.05>138.05 both the PSC and plasma sample are in broad agreement in signal distribution and intensity including the presence of a matrix peak at 2.4mins.

Conclusions

This technology has the potential for a simplified clinical sample collection by the finger prick approach, with future work aimed to evaluate long term sample stability of PSC samples, stored at room temperature. Quantitation of drug metabolites MTIC and AIC also could help provide a measure of sample stability.

References

- Andrasia, M., Bustosb, R., Gaspara, A., Gomez, F.A. & Klekner, A. (2010) **Analysis and stability study of temozolomide using capillary electrophoresis.** *Journal of Chromatography B.* Vol. 878, p1801-1808
- Denny, B.J., Wheelhouse, R.T., Stevens, M.F.G., Tsang, L.L.H., Slack, J.A., (1994) **NMR and molecular modelling investigation of the mechanism of activation of the antitumour drug temozolomide and its Interaction with DNA.** *Biochemistry,* Vol. 33, p9045-9051